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Material: Titan Grad 5

as well as for fabrication of occlusally screwed hybrid-abutment-crowns in the posterior.

**Contraindications:** Primary connections of abutments, Prosthetic angulation to the implant axis of more than 25°, abutment constructions higher than 11 mm measured from implant shoulder, Cast-on technique, Double crowns on implant diameter 3.40 mm, Single tooth restorations with free-end cantilever.

**Processing:**

- Basic procedures for processing are on the authority of the user chairside in practice or labside in the dental laboratory (inLab) taking into account the respective CAD/CAM system.
- The HighLine „C“ Titanium Base is available for all four implant diameters with a gingival height of 0.9 mm.
- For production and further processing of the ceramics part the manufacturers specifications of the CAD/CAM system and of the respective material apply.
- The connection to the implant is guaranteed by a highly precise titanium base with FourByFour® connection. Thus the abutment screw transmits the tear strength to the titanium base during final fixation and not onto the zirconium part of the final abutment.
- Make sure that there is a definite positioning/retention during bonding with the HighLine „C“ titanium base, thus preventing transmission errors.

**IMPORTANT:**

For the scanning process the original Scanbodies of the respective CAD/CAM system have to be used. Under „ADMINISTRATION“ select the implant connection TiBase and the TiBase connection:  
 FX 3.4 for the OT medical implant diameter 3.4 (S)  
 FX 3.8 for the OT medical implant diameter 3.8 (S)  
 FX 4.5 for the OT medical implant diameter 4.1 (L)  
 FX 5.5 for the OT medical implant diameter 5.0 (L)  
 Thereafter please select TiBase in Scanbody type

**Adhesion:** Basically the surfaces of both parts to receive the adhesive should first be blasted carefully with 50µm aluminium oxide. The FourByFour® connection of the titanium base should be protected from damage by screwing (laboratory screw) onto an implant analog prior to blasting and adhesion. The screw entrance of the titanium base should be closed with wax prior to applying adhesive.

For adhesion, please observe the instructions of the adhesive manufacturer. Suitable adhesives are among others:

- Multilink Implant® or Multilink Hybrid Abutment; Ivoclar Vivadent
  - PanaviaTM F 2.0 or PanaviaTM SA Cement; Kuraray
  - RelyX™ Unicem 2 Automix; 3M Espe
- An additional light curing is to be conducted necessarily!

Material: Titanium grade 5

**OT-F<sup>2</sup>/OT-F<sup>3</sup> Scanbody „C“**

**Indications:** The Scanbody „C“ is used as an auxiliary tool for three dimensional registration of the implant in the patients' mouth in order to transfer this in a virtual 3D-model by scanning the correct position. This is the requirement for the design (CAD) and subsequent manufacturing (CAM) of individually machined hybrid abutments.

**Contraindications:** The scanbody „C“ must not be reworked for use as final abutment in the mouth of a patient.

**Processing:** Scanbodies „C“ are color-coded according to the implant diameters and laser marked with “S” or “L” in order to select the compatible Scanbody of the respective CAD/CAM system used in the software.

**IMPORTANT:**

- For intra-oral scan the abutments and the respective screw have to be sterilized prior to use [see enclosed instructions “Preparation (Cleaning, Disinfection and Sterilization) of instruments”, item “impression coping open/closed.”]

**Note:**

The Scanbodies „C“ are manufactured with highest precision for multiple use. Therefore please handle with care and store in a protected way. Also non-obviously recognizable damage or abrasions are leading to faulty scan results. A timely renewal is highly recommended.

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**Introduction:** Illustrations and item numbers referring to the prosthetic components mentioned are contained in detail in the OT-F<sup>2</sup>/OT-F<sup>3</sup> product catalog.

Malocclusion and overloading of the implant abutment connection, the abutments and the complete prosthetic construction should be absolutely avoided.

**Important instructions:**

- All HighLine „C“ titanium bases and Scanbodies „C“ are supplied non-sterile.
- Only new and non-used abutment screws (i.e. not used in the dental laboratory) should be used for final restoration. Titanium bases are supplied with a lab screw mounted and a color-coded final screw. Scanbodies „C“ are supplied with the color coded final screw.
- Implants of  $\varnothing$  3.40/ $\varnothing$  3.80 mm diameter have an internal thread of size M1.6, implants of  $\varnothing$  4.10/ $\varnothing$  5.00 mm an internal thread of M1.8. Please consider this for processing or re-ordering.
- Final screws are color-coded for distinction. Only the final screw may be used for final fixation of the abutments.
- Titanium bases and Scanbodies „C“ have an internal screw thread which prevents a possible loss of the abutment screw during processing.
- Titanium bases have to be fixed with 35 Ncm. Abutment screws may not be cemented or glued within the implant. Scanbodies „C“ should be hand-screwed (15 Ncm).
- While processing of titanium bases the material specific minimum wall thickness has to be observed. The wall thickness may only be reduced to a minimum of 0.5 mm when manufacturing the zirconium part. When using LS2 for crown in the posterior this value is 1,5 mm. Specifications of the manufacturer of the ceramics have to be followed.
- Color markings are provided for Scanbodies „C“ completely, for titanium bases partially in the area of the FourByFour® connection according to the color guide system (green:  $\varnothing$  3.40; yellow:  $\varnothing$  3.80; red:  $\varnothing$  4.10; blue:  $\varnothing$  5.00).
- Connection surfaces of the abutments to the implant may not be blasted or treated. The hazard of an invasion of bacteria as well as weakening of the components must be avoided.

**HighLine „C“ Titanium Base**

**Indications:**

As hybrid abutment for cemented crown and bridge restorations

CE 0482

**OT medical GmbH**

Konsul-Smidt-Straße 8b, 28217 Bremen, Germany  
 Phone +49 421 557161-0, Fax +49 421 557161-95