

Smarter Thinking. Simpler Design.



Genesis

Product Catalog

Dear Customer,

Thank you for allowing Keystone Dental Inc. the opportunity to provide solutions which we hope meet or exceed your expectations.

Keystone Dental Inc. produces innovative, high quality dental implants, accessories, and biologics and distributes these products globally. Our headquarters are in Burlington, Massachusetts and we have state-of-the-art manufacturing facilities in both Burlington and Irvine, California.

Our range of implant products is broad, fulfilling both niche and standard requirements for all specialties. If you wish, at any time, to discuss Keystone Dental products, ongoing developments, and future research, or if you have any suggestions which you feel would help us serve you better, please call me directly at (781) 328-3387 or Toll Free at (866) 902-9272.

Again, thank you for allowing Keystone Dental, Inc. to satisfy your implant requirements, and I look forward to speaking with you in the future.

A handwritten signature in white ink that reads "Michael Kehoe". The signature is written in a cursive, flowing style.

Michael A. Kehoe
Keystone Dental, Inc.
President & CEO




Genesis

Ø3.5/3.8 mm Genesis Implants	4
Ø3.5/3.8 mm Prosthetic Options	5
Ø3.5/3.8 mm Prosthetic Components	6
Ø4.5 mm Genesis Implants	10
Ø4.5 mm Prosthetic Options	11
Ø4.5 mm Prosthetic Components	12
Ø5.5/6.5 mm Genesis Implants	16
Ø5.5/6.5 mm Prosthetic Options	17
Ø5.5/6.5 mm Prosthetic Components	18
LOCATOR® Prosthetic Options	23
LOCATOR R-Tx™ Abutment	24
LOCATOR® Abutment	25
Genesis Surgical Kit	26
Surgical Motor	30
Terms & Conditions	31

Product specifications are subject to change without notice.
Items illustrated are not to scale.



Genesis Implants



Catalog #

G21050
G21051
G21052
G21053

G21130
G21131
G21132
G21133
G21135

Implant Size

Ø3.5 x 10 mm Length
Ø3.5 x 11.5 mm Length
Ø3.5 x 13 mm Length
Ø3.5 x 16 mm Length

Ø3.8 x 8.5 mm Length
Ø3.8 x 10 mm Length
Ø3.8 x 11.5 mm Length
Ø3.8 x 13 mm Length
Ø3.8 x 16 mm Length

Cover Screw



Catalog #

G21100

Description

Ø3.5/3.8 mm

Healing Abutments



Catalog #

G21103
G21104
G21105
G21106

Description

Ø3.5/3.8 mm
Ø3.5/3.8 mm
Ø3.5/3.8 mm
Ø3.5/3.8 mm

Flare


4.0 mm
4.0 mm
5.0 mm
5.0 mm

Cuff Height

3.0 mm
5.0 mm
3.0 mm
5.0 mm









The Ø3.5 & Ø3.8 mm Genesis can be restored with Ø4.0 and Ø5.0 mm flared prosthetic components. An Abutment Screw is included with all TiLobe abutments.

Ø3.5 mm Implant available in 2016.






Prosthetic Components
Internal TiLobe

	Catalog #	Description	Flare	
	G71100	Ø3.5/3.8 mm	4.0 mm	
	G71101	Ø3.5/3.8 mm	5.0 mm	
	45904K	Ø3.5/3.8 mm	4.0 mm	
	45906K	Ø3.5/3.8 mm	5.0 mm	
	45905K	Ø3.5/3.8 mm	4.0 mm	
	45907K	Ø3.5/3.8 mm	5.0 mm	
	G71107	Ø3.5/3.8 mm, 1-Pack		
	G71108	Ø3.5/3.8 mm, 10-Pack		
	G71130	Ø3.5/3.8 mm, Locking	4.0 mm	
	G71131	Ø3.5/3.8 mm, Non-Locking	4.0 mm	
	G71132	Ø3.5/3.8 mm, Locking	5.0 mm	
	G71133	Ø3.5/3.8 mm, Non-Locking	5.0 mm	
	G71144	Ø3.5/3.8 mm, Locking		
	G71145	Ø3.5/3.8 mm, Non-Locking		
 Cuff Height (Buccal)	G71200	Ø3.5/3.8 mm, Straight	4.0 mm	1.0 mm
	G71201	Ø3.5/3.8 mm, Straight	4.0 mm	2.0 mm
	G71202	Ø3.5/3.8 mm, 15° Angled	4.0 mm	1.0 mm
	G71203	Ø3.5/3.8 mm, 15° Angled	4.0 mm	2.0 mm
	G71204	Ø3.5/3.8 mm, Straight	5.0 mm	1.0 mm
	G71205	Ø3.5/3.8 mm, Straight	5.0 mm	3.0 mm
	G71206	Ø3.8 mm, 15° Angled	5.0 mm	1.0 mm
	G71207	Ø3.8 mm, 15° Angled	5.0 mm	3.0 mm
	G71231	Gold/Plastic Sleeve, Locking		1.0 mm
	G71232	Gold/Plastic Sleeve, Non-Locking		1.0 mm



	Replacement Final Abutment Screw	<i>Catalog #</i>	<i>Description</i>		
		G71125	1-Pack		
	Laboratory Screw	<i>Catalog #</i>	<i>Description</i>		
		G71127	10-Pack		
	Quick-Packs	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
		G71400	Ø3.5/3.8 mm	5.0 mm	1.0 mm
		G71401	Ø3.5/3.8 mm	5.0 mm	2.0 mm
		G71402	Ø3.5/3.8 mm	5.0 mm	3.0 mm
		G71403	Ø3.5/3.8 mm	5.0 mm	4.0 mm
	Quick-Abutment	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
		G71434	Ø3.5/3.8 mm	5.0 mm	1.0 mm
		G71435	Ø3.5/3.8 mm	5.0 mm	2.0 mm
		G71436	Ø3.5/3.8 mm	5.0 mm	3.0 mm
		G71437	Ø3.5/3.8 mm	5.0 mm	4.0 mm
	Quick-Abutment Quick-Cap®	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	
		G71428	Ø3.5/3.8 mm, 1-Pack	5.0 mm	
	Quick-Abutment Temporary/Healing Cap	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	
		G71462	Ø3.5/3.8 mm	5.0 mm	
	Quick-Abutment Analog	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	
		G71468	Ø3.5/3.8 mm	5.0 mm	
	Quick-Abutment Plastic Sleeve	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	
		G71478	Ø3.5/3.8 mm, Locking	5.0 mm	
		G71479	Ø3.5/3.8 mm, Non-Locking	5.0 mm	
	Quick-Abutment Ceramic Coping	<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	
		G71475	Ø3.5/3.8 mm	5.0 mm	




Ti-Base Abutment (Short)

Catalog #

Description

45806K

Ø3.5/3.8 mm




Pink Ti-Base Abutment (Short)

Catalog #

Description

45007K

Ø3.5/3.8 mm




Ti-Base Abutment (Long)

Catalog #

Description

45901K

Ø3.5/3.8 mm



Multi Unit Abutment

Catalog #

Description

Cuff Height

45312K

Ø3.5/3.8 mm, Straight

1.0 mm

45313K

Ø3.5/3.8 mm, Straight

2.0 mm

45314K

Ø3.5/3.8 mm, Straight

3.0 mm

45321K

Ø3.5/3.8 mm, Angled, 17°

3.0 mm

45322K

Ø3.5/3.8 mm, Angled, 17°

4.0 mm

Multi Unit Abutment Healing Cap

Catalog #

Description



45330K

Ø3.5/3.8 mm



Multi Unit Abutment Impression Post

Catalog #

Description

45327K

Open Tray, Ø3.5/3.8 mm with screw

45328K

Closed Tray, Ø3.5/3.8 mm

Multi Unit Abutment Analog

Catalog #

Description



45329K

Ø3.5/3.8 mm



**Multi Unit Abutment
Ti Temporary Sleeve**

Catalog #

Description

45331K

Ø3.5/3.8 mm



**Multi Unit Abutment
Plastic Sleeve**

Catalog #

Description

45332K

Ø3.5/3.8 mm



**Multi Unit Abutment
Prosthetic Screw**

Catalog #

Description

45333K

1-Pack




**Replacement Abutment
Screw for Multi Unit
Angled Abutment**

Catalog #

Description

45334K

1-Pack



Genesis Implants

Catalog #

- G21137**
- G21138**
- G21139**
- G21140**
- G21142**

Implant Size

- Ø4.5 x 8.5 mm Length
- Ø4.5 x 10 mm Length
- Ø4.5 x 11.5 mm Length
- Ø4.5 x 13 mm Length
- Ø4.5 x 16 mm Length



Cover Screw

Catalog #

- G21101**

Description

- Ø4.5 mm



Healing Abutments

Catalog #

- G21107**
- G21108**
- G21109**
- G21110**

Description

- Ø4.5 mm
- Ø4.5 mm
- Ø4.5 mm
- Ø4.5 mm


Flare

- 5.0 mm
- 5.0 mm
- 6.0 mm
- 6.0 mm

Cuff Height

- 3.0 mm
- 5.0 mm
- 3.0 mm
- 5.0 mm

The Ø4.5 mm Genesis can be restored with Ø5.0 and Ø6.0 mm flared prosthetic components.
An Abutment Screw is included with all TiLobe abutments.





Impression Post

Catalog #

Description

Flare

G71102
G71103

Ø4.5 mm
Ø4.5 mm

5.0 mm
6.0 mm



**Open Tray
Impression Post**

Catalog #

Description

Flare

45908K
45910K

Ø4.1/4.5 mm
Ø4.1/4.5 mm

5.0 mm
6.0 mm



**Closed Tray
Impression Post**

Catalog #

Description

Flare

45909K
45911K

Ø4.1/4.5 mm
Ø4.1/4.5 mm

5.0 mm
6.0 mm



Implant Analog

Catalog #

Description

G71109
G71110

Ø4.5 mm, 1-Pack
Ø4.5 mm, 10-Pack



Temporary Abutment

Catalog #


Description

Flare

G71134
G71135
G71136
G71137

Ø4.5 mm, Locking
Ø4.5 mm, Non-Locking
Ø4.5 mm, Locking
Ø4.5 mm, Non-Locking

5.0 mm
5.0 mm
6.0 mm
6.0 mm



Ti Temporary Abutment

Catalog #

Description

G71146
G71147

Ø4.5 mm, Locking
Ø4.5 mm, Non-Locking



Cuff Height
(Buccal)

**Esthetic Contour Ti
Abutment**

Catalog #

Description

Flare

Cuff Height

G71208
G71209
G71210
G71211

Ø4.5 mm, Straight
Ø4.5 mm, Straight
Ø4.5 mm, 15° Angled
Ø4.5 mm, 15° Angled

5.0 mm
5.0 mm
5.0 mm
5.0 mm

0.5 mm
1.5 mm
0.5 mm
1.5 mm

G71212
G71213
G71214
G71215

Ø4.5 mm, Straight
Ø4.5 mm, Straight
Ø4.5 mm, 15° Angled
Ø4.5 mm, 15° Angled

6.0 mm
6.0 mm
6.0 mm
6.0 mm

1.0 mm
3.0 mm
1.0 mm
3.0 mm



UCLA Abutment

Catalog #

Description

Cuff Height

G71233
G71234

Gold/Plastic Sleeve, Locking
Gold/Plastic Sleeve, Non-Locking

1.0 mm
1.0 mm



**Replacement
Final Abutment Screw**

Catalog #

Description

G71125

1-Pack



Laboratory Screw

Catalog #

Description

G71127

10-Pack

Quick-Packs

Catalog #

Description

Flare

Cuff Height



G71404

Ø4.5 mm

5.0 mm

0.5 mm

G71405

Ø4.5 mm

5.0 mm

1.5 mm

G71406

Ø4.5 mm

5.0 mm

2.5 mm

G71407

Ø4.5 mm

5.0 mm

4.0 mm

G71408

Ø4.5 mm

6.0 mm

0.5 mm

G71409

Ø4.5 mm

6.0 mm

1.5 mm

G71410

Ø4.5 mm

6.0 mm

2.5 mm

G71411

Ø4.5 mm

6.0 mm

4.0 mm


Quick-Abutment

Catalog #

Description

Flare

Cuff Height



G71438

Ø4.5 mm

5.0 mm

0.5 mm

G71439

Ø4.5 mm

5.0 mm

1.5 mm

G71440

Ø4.5 mm

5.0 mm

2.5 mm

G71441

Ø4.5 mm

5.0 mm

4.0 mm

G71442

Ø4.5 mm

6.0 mm

0.5 mm

G71443

Ø4.5 mm

6.0 mm

1.5 mm

G71444

Ø4.5 mm

6.0 mm

2.5 mm

G71445

Ø4.5 mm

6.0 mm

4.0 mm

**Quick-Abutment
Quick-Cap®**

Catalog #

Description

Flare



G71429

Ø4.5 mm, 1-Pack

5.0 mm

G71430

Ø4.5 mm, 1-Pack

6.0 mm

**Quick-Abutment
Temporary/Healing Cap**

Catalog #

Description

Flare



G71463

Ø4.5 mm

5.0 mm

G71464

Ø4.5 mm

6.0 mm

**Quick-Abutment
Analog**

Catalog #

Description

Flare



G71469

Ø4.5 mm

5.0 mm

G71470

Ø4.5 mm

6.0 mm



**Quick-Abutment
Plastic Sleeve**

Catalog #

Description

Flare

G71480

Ø4.5 mm, Locking

5.0 mm

G71481

Ø4.5 mm, Non-Locking

5.0 mm

G71482

Ø4.5 mm, Locking

6.0 mm

G71483

Ø4.5 mm, Non-Locking

6.0 mm



**Quick-Abutment
Ceramic Coping**

Catalog #

Description

Flare

G71476

Ø4.5 mm

5.0 mm



Multi Unit Abutment

Catalog #

Description

Cuff Height

45315K

Ø4.1/4.5 mm, Straight

1.0 mm

45316K

Ø4.1/4.5 mm, Straight

2.0 mm

45317K

Ø4.1/4.5 mm, Straight

3.0 mm

45318K

Ø4.1/4.5 mm, Straight

4.0 mm

45323K

Ø4.1/4.5 mm, Angled, 17°

3.0 mm

45324K

Ø4.1/4.5 mm, Angled, 17°

4.0 mm

45325K

Ø4.1/4.5 mm, Angled, 30°

4.0 mm

45326K

Ø4.1/4.5 mm, Angled, 30°

5.0 mm



**Multi Unit Abutment
Healing Cap**

Catalog #

Description

45330K

Ø4.1/4.5 mm



**Multi Unit Abutment
Impression Post**

Catalog #

Description

45327K

Open Tray, Ø4.1/4.5 mm with screw

45328K

Closed Tray, Ø4.1/4.5 mm



**Multi Unit Abutment
Analog**

Catalog #

Description

45329K

Ø4.1/4.5 mm



**Multi Unit Abutment
Ti Temporary Sleeve**

Catalog #

Description

45331K

Ø4.1/4.5 mm



**Multi Unit Abutment
Plastic Sleeve**

Catalog #

Description

45332K

Ø4.1/4.5 mm



**Multi Unit Abutment
Prosthetic Screw**

Catalog #

Description

45333K

1-Pack



**Replacement Abutment
Screw for Multi Unit
Angled Abutment**

Catalog #

Description

45334K

1-Pack




Ti-Base Abutment (Short)

Catalog #

Description

45807K

Ø4.1/4.5 mm




**Pink Ti-Base
Abutment (Short)**

Catalog #

Description

45008K

Ø4.1/4.5 mm




Ti-Base Abutment (Long)

Catalog #

Description

45902K

Ø4.1/4.5 mm



Genesis Implants

Catalog #

Implant Size

- G21144** Ø5.5 x 8.5 mm Length
- G21145** Ø5.5 x 10 mm Length
- G21146** Ø5.5 x 11.5 mm Length
- G21147** Ø5.5 x 13 mm Length

Catalog #

Implant Size

- G21150** Ø6.5 x 8.5 mm Length
- G21151** Ø6.5 x 10 mm Length
- G21152** Ø6.5 x 11.5 mm Length
- G21153** Ø6.5 x 13 mm Length

Cover Screw



Catalog #

Description

- G21102** Ø5.5/6.5 mm

Healing Abutments



Catalog #


Description

Flare

Cuff Height

- | | | | |
|---------------|-------------|--------|--------|
| G21111 | Ø5.5 mm | 5.0 mm | 3.0 mm |
| G21112 | Ø5.5 mm | 5.0 mm | 5.0 mm |
| G21113 | Ø5.5/6.5 mm | 6.0 mm | 3.0 mm |
| G21114 | Ø5.5/6.5 mm | 6.0 mm | 5.0 mm |
| G21296 | Ø5.5/6.5 mm | 7.0 mm | 3.0 mm |
| G21297 | Ø5.5/6.5 mm | 7.0 mm | 5.0 mm |
| G21298 | Ø5.5/6.5 mm | 8.0 mm | 3.0 mm |
| G21299 | Ø5.5/6.5 mm | 8.0 mm | 5.0 mm |

The Ø5.5 and Ø6.5 mm Genesis can be restored with Ø5.0 and Ø6.0 mm flared prosthetic components. An Abutment Screw is included with all TiLobe abutments.





Impression Post

Catalog #

Description

Flare

G71104
G71105

Ø5.5 mm
Ø5.5/6.5 mm

5.0 mm
6.0 mm



Impression Post

Catalog #

Description

Flare

45916K
45918K

Ø5.5/6.5 mm
Ø5.5/6.5 mm

7.0 mm
8.0 mm



Impression Post

Catalog #

Description

Flare

45917K
45919K

Ø5.5/6.5 mm
Ø5.5/6.5 mm

7.0 mm
8.0 mm



Implant Analog

Catalog #

Description

G71111
G71112
G71113
G71114

Ø5.5 mm, 1-Pack
Ø5.5 mm, 10-Pack
Ø6.5 mm 1-Pack
Ø6.5 mm 10-Pack



Temporary Abutment

Catalog #


Description

Flare

G71138
G71139
G71140
G71141

Ø5.5 mm, Locking
Ø5.5 mm, Non-Locking
Ø5.5/6.5 mm, Locking
Ø5.5/6.5 mm, Non-Locking

5.0 mm
5.0 mm
6.0 mm
6.0 mm



Ti Temporary Abutment

Catalog #

Description

G71148
G71149

Ø5.5/6.5 mm, Locking
Ø5.5/6.5 mm, Non-Locking



Cuff Height
(Buccal)

Esthetic Contour Ti Abutment

Catalog #

Description

Flare

Cuff Height

G71216
G71217
G71218
G71219

Ø5.5 mm, Straight
Ø5.5 mm, Straight
Ø5.5 mm, 15° Angled
Ø5.5 mm, 15° Angled

5.0 mm 0.5 mm
5.0 mm 1.5 mm
5.0 mm 0.5 mm
5.0 mm 1.5 mm

G71220
G71221
G71222

Ø5.5/6.5 mm, Straight
Ø5.5/6.5 mm, Straight
Ø5.5/6.5 mm, 15° Angled

6.0 mm 1.0 mm
6.0 mm 3.0 mm
6.0 mm 1.0 mm



UCLA Abutment

Catalog #

Description

Cuff Height

G71235
G71236

Gold/Plastic Sleeve, Locking
Gold/Plastic Sleeve, Non-Locking

1.0 mm
1.0 mm

**Replacement
Final Abutment Screw**



<i>Catalog #</i>	<i>Description</i>
G71125	1-Pack

Laboratory Screw




<i>Catalog #</i>	<i>Description</i>
G71127	10-Pack

Quick-Packs



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
G71412	Ø5.5 mm	5.0 mm	0.5 mm
G71413	Ø5.5 mm	5.0 mm	1.5 mm
G71414	Ø5.5 mm	5.0 mm	2.5 mm
G71415	Ø5.5 mm	5.0 mm	4.0 mm
G71416	Ø5.5/6.5 mm	6.0 mm	0.5 mm
G71417	Ø5.5/6.5 mm	6.0 mm	1.5 mm
G71418	Ø5.5/6.5 mm	6.0 mm	2.5 mm
G71419	Ø5.5/6.5 mm	6.0 mm	4.0 mm

Quick-Abutment



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>	<i>Cuff Height</i>
G71446	Ø5.5 mm	5.0 mm	0.5 mm
G71447	Ø5.5 mm	5.0 mm	1.5 mm
G71448	Ø5.5 mm	5.0 mm	2.5 mm
G71449	Ø5.5 mm	5.0 mm	4.0 mm
G71450	Ø5.5/6.5 mm	6.0 mm	0.5 mm
G71451	Ø5.5/6.5 mm	6.0 mm	1.5 mm
G71452	Ø5.5/6.5 mm	6.0 mm	2.5 mm
G71453	Ø5.5/6.5 mm	6.0 mm	4.0 mm

**Quick-Abutment
Quick-Cap®**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
G71431	Ø5.5 mm, 1-Pack	5.0 mm
G71432	Ø5.5/6.5 mm, 1-Pack	6.0 mm

**Quick-Abutment
Temporary/Healing Cap**



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
G71465	Ø5.5 mm	5.0 mm
G71466	Ø5.5/6.5 mm	6.0 mm

Quick-Abutment Analog



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
G71471	Ø5.5 mm	5.0 mm
G71472	Ø5.5/6.5 mm	6.0 mm



Quick-Abutment Plastic Sleeve




<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
G71484	Ø5.5 mm, Locking	5.0 mm
G71485	Ø5.5 mm, Non-Locking	5.0 mm
G71486	Ø5.5/6.5 mm, Locking	6.0 mm
G71488	Ø5.5/6.5 mm, Non-Locking	6.0 mm

Quick-Abutment Ceramic Coping



<i>Catalog #</i>	<i>Description</i>	<i>Flare</i>
G71477	Ø5.5 mm	6.0 mm

Multi Unit Abutment



<i>Catalog #</i>	<i>Description</i>	<i>Cuff Height</i>
45319K	Ø5.0/5.5/6.5 mm, Straight	1.0 mm
45320K	Ø5.0/5.5/6.5 mm, Straight	2.0 mm

Multi Unit Abutment Healing Cap



<i>Catalog #</i>	<i>Description</i>
45330K	Ø5.0/5.5/6.5 mm

Multi Unit Abutment Impression Post



<i>Catalog #</i>	<i>Description</i>
45327K	Open Tray, Ø5.0/5.5/6.5 mm with screw
45328K	Closed Tray, Ø5.0/5.5/6.5 mm

Multi Unit Abutment Analog



<i>Catalog #</i>	<i>Description</i>
45329K	Ø5.0/5.5/6.5 mm

Multi Unit Abutment Ti Temporary Sleeve



<i>Catalog #</i>	<i>Description</i>
45331K	Ø5.0/5.5/6.5 mm



**Multi Unit Abutment
Plastic Sleeve**

Catalog #

Description

45332K

Ø5.0/5.5/6.5 mm



**Multi Unit Abutment
Prosthetic Screw**

Catalog #

Description

45333K

1-Pack




**Replacement Abutment
Screw for Multi Unit
Angled Abutment**

Catalog #

Description

45334K

1-Pack




Ti-Base Abutment (Short)

Catalog #

Description

45808K

Ø5.0/5.5/6.5 mm



**Pink Ti-Base
Abutment (Short)**

Catalog #

Description

45009K

Ø5.0/5.5/6.5 mm




Ti-Base Abutment (Long)

Catalog #

Description

45903K

Ø5.0/5.5/6.5 mm



LOCATOR R-Tx™ Abutment

Internal TiLobe

LOCATOR R-Tx™ Abutment



Catalog

LOCR-TX 31200-01
LOCR-TX 31200-02
LOCR-TX 31200-03
LOCR-TX 31200-04
LOCR-TX 31200-05
LOCR-TX 31200-06

Description

Ø3.5/3.8 mm
 Ø3.5/3.8 mm
 Ø3.5/3.8 mm
 Ø3.5/3.8 mm
 Ø3.5/3.8 mm
 Ø3.5/3.8 mm

Cuff Height

1.0 mm
 2.0 mm
 3.0 mm
 4.0 mm
 5.0 mm
 6.0 mm

LOCR-TX 31201-01
LOCR-TX 31201-02
LOCR-TX 31201-03
LOCR-TX 31201-04
LOCR-TX 31201-05
LOCR-TX 31201-06

Ø4.1/4.5 mm
 Ø4.1/4.5 mm
 Ø4.1/4.5 mm
 Ø4.1/4.5 mm
 Ø4.1/4.5 mm
 Ø4.1/4.5 mm

1.0 mm
 2.0 mm
 3.0 mm
 4.0 mm
 5.0 mm
 6.0 mm

LOCR-TX 31202-01
LOCR-TX 31202-02
LOCR-TX 31202-03
LOCR-TX 31202-04
LOCR-TX 31202-05
LOCR-TX 31202-06

Ø5.0/5.5/6.5 mm
 Ø5.0/5.5/6.5 mm
 Ø5.0/5.5/6.5 mm
 Ø5.0/5.5/6.5 mm
 Ø5.0/5.5/6.5 mm
 Ø5.0/5.5/6.5 mm

1.0 mm
 2.0 mm
 3.0 mm
 4.0 mm
 5.0 mm
 6.0 mm

LOCATOR R-Tx™ Abutment Accessories



Catalog

LOCR-TX 30013-01
LOCR-TX 30001-01
LOCR-TX 30002-01
LOCR-TX 30003-01
LOCR-TX 30004-01
LOCR-TX 30005-01
LOCR-TX 30014-01
LOCR-TX 30015-01
LOCR-TX 30016-01
LOCR-TX 30017-01
LOCR-TX 30021-01
LOCR-TX 04936
LOCR-TX 08007
LOCR-TX 08008
LOCR-TX 30012-01
LOCR-TX 30018-01

Description

Denture Attachment Processing Assembly, 4-Pack
 Zero Retention Insert, 4-Pack
 Low Retention Insert, 4-Pack
 Medium Retention Insert, 4-Pack
 High Retention Insert, 4-Pack
 Extra-High Retention Insert, 4-Pack
 3.35 mm Abutment Analog, 4-Pack
 4.0 mm Abutment Analog, 4-Pack
 5.0 mm Abutment Analog, 4-Pack
 Impression Coping, 4-Pack
 Retention Insert Tool
 .050"/1.25 mm Hex, Torque Wrench
 .050"/1.25 mm Hex Driver with Thumb Knob
 .050"/1.25 mm Hex, Latch Lock Driver
 Processing Insert, 4-Pack
 Processing Spacer, 4-Pack



LOCATOR® Abutment Assembly




<i>Catalog #</i>	<i>Description</i>	<i>Cuff Height</i>
G71500	Ø3.5/3.8 mm	1.0 mm
G71501	Ø3.5/3.8 mm	2.0 mm
G71502	Ø3.5/3.8 mm	3.0 mm
G71503	Ø3.5/3.8 mm	4.0 mm
G71504	Ø3.5/3.8 mm	5.0 mm
G71505	Ø3.5/3.8 mm	6.0 mm
G71506	Ø4.5 mm	1.0 mm
G71507	Ø4.5 mm	2.0 mm
G71508	Ø4.5 mm	3.0 mm
G71509	Ø4.5 mm	4.0 mm
G71510	Ø4.5 mm	5.0 mm
G71511	Ø4.5 mm	6.0 mm
G71512	Ø5.5 mm	1.0 mm
G71513	Ø5.5 mm	2.0 mm
G71514	Ø5.5 mm	3.0 mm
G71515	Ø5.5 mm	4.0 mm
G71516	Ø5.5 mm	5.0 mm
G71517	Ø5.5 mm	6.0 mm


LOCATOR® Abutment Accessories



<i>Catalog #</i>	<i>Description</i>
L8519-2	Denture Cap with Black Processing Male, 2-Pack
G71518	LOCATOR® Replacement Male, White, Standard Retention, 4-Pack
G71519	LOCATOR® Replacement Male, Pink, Light Retention, 4-Pack
G71520	LOCATOR® Replacement Male, Blue, (ELR) Extra Light Retention, 4-Pack
G71521	LOCATOR® Replacement Male, Green, (ER) Extended Range, 4-Pack
G71522	LOCATOR® Female Analog, 4-Pack
G71523	LOCATOR® Impression Coping, 4-Pack
G71524	LOCATOR® Core Tool
G71525	LOCATOR® Torque Wrench Driver




Genesis
Surgical Kit



Catalog #	Surgical Options
G31000	Surgical Cassette (empty)
G31001	A. Genesis Surgical Kit for Ø3.5-Ø4.5 mm
G31002	B. Genesis Drill Pack for Ø5.5/Ø6.5 mm

Catalog #	Description	Catalog #	Description
G31001	A. Genesis Surgical Kit for Ø3.5-Ø4.5 mm <i>Includes:</i>	G31036	Implant Depth Gauge 4.5 x 8.5-10 mm
G31000	Surgical Cassette	G31037	Implant Depth Gauge 4.5 x 11.5-13 mm
G21244	2.0 mm Full Depth Twist Drill	G31038	Implant Depth Gauge 4.5 x 16 mm
G21211	Parallel Pin 4-Pack	G21290	3.5/3.8 Implant Driver, Latch, Short
G31008	Final Tapered Drill 3.5 x 8.5 mm	G21291	3.5/3.8 Implant Driver, Latch, Long
G31009	Final Tapered Drill 3.5 x 10 mm	G21292	4.5 Implant Driver, Latch, Short
G31010	Final Tapered Drill 3.5 x 11.5 mm	G21293	4.5 Implant Driver, Latch, Long
G31011	Final Tapered Drill 3.5 x 13 mm	G21274	3.5/3.8 Implant Driver, Ratchet
G31012	Final Tapered Drill 3.5 x 16 mm	G21275	4.5 Implant Driver, Ratchet
G31013	Final Tapered Drill 3.8 x 8.5 mm	K15695	Surgical Ratchet/Torque Wrench
G31014	Final Tapered Drill 3.8 x 10 mm	15224K	Ratchet Adapter
G31015	Final Tapered Drill 3.8 x 11.5 mm	G71612	Quad Driver, Torque Wrench, Short
G31016	Final Tapered Drill 3.8 x 13 mm	G71613	Quad Driver, Torque Wrench, Long
G31017	Final Tapered Drill 3.8 x 16 mm	G71608	Quad Driver, Swivel Head, Short
G31018	Final Tapered Drill 4.5 x 8.5 mm	G71609	Quad Driver, Swivel Head, Long
G31019	Final Tapered Drill 4.5 x 10 mm	G21204	Drill Extender
G31020	Final Tapered Drill 4.5 x 11.5 mm	G21277	Conversion Handle, Latch Type
G31021	Final Tapered Drill 4.5 x 13 mm	G21278	Conversion Handle, Square Type
G31022	Final Tapered Drill 4.5 x 16 mm	100001-EN01	Surgical Manual
G31031	Implant Depth Gauge 3.5 x 10-11.5 mm	100009-EN01	Drill Usage Chart
G31032	Implant Depth Gauge 3.5 x 13-16 mm	100005-EN01	Wall Chart
G31033	Implant Depth Gauge 3.8 x 8.5-10 mm	100027-EN01	Implant Driver, Latch IFU
G31034	Implant Depth Gauge 3.8 x 11.5 -13 mm	MK40091	Genesis Folder
G31035	Implant Depth Gauge 3.8 x 16 mm		



<i>Catalog #</i>	<i>Description</i>
G31002	B. Genesis Drill Pack for Ø5.5/6.5 mm <i>Includes:</i>
G31023	Final Tapered Drill 5.5 x 8.5 mm
G31024	Final Tapered Drill 5.5 x 10 mm
G31025	Final Tapered Drill 5.5 x 11.5 mm
G31026	Final Tapered Drill 5.5 x 13 mm
G31027	Final Tapered Drill 6.5 x 8.5 mm
G31028	Final Tapered Drill 6.5 x 10 mm
G31029	Final Tapered Drill 6.5 x 11.5 mm
G31030	Final Tapered Drill 6.5 x 13 mm
G31039	Implant Depth Gauge 5.5 x 8.5-10 mm
G31040	Implant Depth Gauge 5.5 x 11.5-13 mm
G31041	Implant Depth Gauge 6.5 x 8.5-10 mm
G31042	Implant Depth Gauge 6.5 x 11.5-13 mm
G21294	5.5/6.5 Implant Driver, Latch, Short
G21295	5.5/6.5 Implant Driver, Latch, Long
G21276	5.5/6.5 Implant Driver, Ratchet
100027-EN01	Implant Driver, Latch IFU

Initial Drill




<i>Catalog #</i>	<i>Description</i>
G21243	Spade Drill
G21244	2.0 mm Full Depth Twist Drill

Final Tapered Drills



<i>Catalog #</i>	<i>Description</i>
G31008	Ø3.5 x 8.5 mm
G31009	Ø3.5 x 10 mm
G31010	Ø3.5 x 11.5 mm
G31011	Ø3.5 x 13 mm
G31012	Ø3.5 x 16 mm
G31013	Ø3.8 x 8.5 mm
G31014	Ø3.8 x 10 mm
G31015	Ø3.8 x 11.5 mm
G31016	Ø3.8 x 13 mm
G31017	Ø3.8 x 16 mm
G31018	Ø4.5 x 8.5 mm
G31019	Ø4.5 x 10 mm
G31020	Ø4.5 x 11.5 mm
G31021	Ø4.5 x 13 mm
G31022	Ø4.5 x 16 mm
G31023	Ø5.5 x 8.5 mm
G31024	Ø5.5 x 10 mm
G31025	Ø5.5 x 11.5 mm
G31026	Ø5.5 x 13 mm



Final Tapered Drills (cont.)



Catalog #

- G31027**
- G31028**
- G31029**
- G31030**

Description

- Ø6.5 x 8.5 mm
- Ø6.5 x 10 mm
- Ø6.5 x 11.5 mm
- Ø6.5 x 13 mm

Depth Gauges



Catalog #

- G31031**
- G31032**

Description

- Ø3.5 x 10-11.5 mm
- Ø3.5 x 13-16 mm

- G31033**
- G31034**
- G31035**

- Ø3.8 x 8.5-10 mm
- Ø3.8 x 11.5-13 mm
- Ø3.8 x 16 mm

- G31036**
- G31037**
- G31038**

- Ø4.5 x 8.5-10 mm
- Ø4.5 x 11.5-13 mm
- Ø4.5 x 16 mm

- G31039**
- G31040**

- Ø5.5 x 8.5-10 mm
- Ø5.5 x 11.5-13 mm

- G31041**
- G31042**

- Ø6.5 x 8.5-10 mm
- Ø6.5 x 11.5-13 mm

Surgical Taps



Catalog #

- G31045**
- G31046**

Description

- Ø3.5 x 10-11.5 mm
- Ø3.5 x 13-15 mm

- G21008**

- Ø3.8 x 11.5-16 mm

- G21009**
- G21010**

- Ø4.5 x 10-13 mm
- Ø4.5 x 16 mm

- G21011**
- G21012**

- Ø5.5 x 8.5 mm
- Ø5.5 x 10-13 mm

- G21013**
- G21014**

- Ø6.5 x 8.5 mm
- Ø6.5 x 10-13 mm

Parallel Pin



Catalog #

- G21211**

Description


- 4-Pack

Drill Extender



Catalog #

- G21204**



Implant Driver, Latch



Catalog #

G21290
G21291

Description

Ø3.5/3.8 mm, Short
Ø3.5/3.8 mm, Long


G21292
G21293

Ø4.5 mm, Short
Ø4.5 mm, Long

G21294
G21295

Ø5.5/6.5 mm, Short
Ø5.5/6.5 mm, Long

Implant Driver, Ratchet



Catalog #

G21274
G21275
G21276

Description

Ø3.5/3.8 mm
Ø4.5 mm
Ø5.5/6.5 mm

**Surgical Ratchet/
Torque Wrench**



Catalog #

K15695

Description

Surgical Ratchet/Torque Wrench

Ratchet Adapter



Catalog #

15224K

Conversion Handles



Catalog #

G21277
G21278

Description

Conversion Handle, Latch Type
Conversion Handle, Square Type

Quad Drivers



Catalog #

G71608
G71609

Description

Swivel Head, Short
Swivel Head, Long



G71610
G71611

Latch Type, Short
Latch Type, Long



G71612
G71613

Torque Wrench, Short
Torque Wrench, Long



Implantmed Surgical Motor & Accessories



<i>Catalog #</i>	<i>Description</i>
WH-16929001	SI-915 Implantmed Surgical Motor Unit
WH-10207550	WI-75 E/KM Contra-Angle 20:1 Handpiece for implant surgery (15-2,000 rpm)
WH-30065000	WS-75 Contra-Angle 20:1 Dismantlable Handpiece for implant surgery (15-2,000 rpm)
WH-30032000	WS-75 L G Contra-Angle 20:1 Dismantlable Handpiece with Mini LED+ for implant surgery (15-2,000 rpm)
WH-30056000	S-11 1:1 High Speed Straight Handpiece for oral surgery (300-40,000 rpm)
WH-04363600	Irrigation Tubing and Spike Assembly, 6-Pack, Disposable (for use with Implantmed series & Elcomed SA-310 motors)
WH-00929300	Irrigation Tubing with Luer Lock, Green (10-Pack)
WH-04019000	Irrigation Tubing Clips, 5-Pack (for use with Implantmed series motors to hold tubing to motor cable)
WH-02610500	Irrigation Needle/Canula (for use with internally irrigated drills)
WH-06016400	Y-Connector Complete
WH-06946300	Irrigation Spray Clip (new series), Left Side, E/KM, Green (3-Pack)
WH-06948300	Irrigation Spray Clip (new series), Right Side, E/KM, Green (3-Pack)
WH-06948400	Irrigation Spray Clip (new series), Left Side, Green (3-Pack)
WH-06949500	Irrigation Spray Clip (new series), Right Side, Green (3-Pack)
WH-04757100	Irrigation Spray Clip (old series), E/KM, Black (3-Pack)
WH-05051000	Irrigation Spray Clip (old series), Black (3-Pack)
WH-02139800	Bur Testing Gauge (to check all 2.35 mm & 1.6 mm burs)
WH-10940011	W&H F-1 Service Oil, 400 ml (Special shipping & handling required [aerosol can]. Contact Customer Service for shipping rates.)
WH-02038200	Spray Cap for Service Oil (reusable replacement)
WH-04013500	Sterilization Case

May not be available in all markets.

Terms and Conditions

TERMS AND CONDITIONS OF SALE. The following are the terms and conditions by which Keystone Dental, Inc. sells its products. As used herein, the term “Keystone” means Keystone Dental, Inc. and its applicable affiliates; the term “Products” means the Keystone products and services described in this invoice, and “Purchaser” means the person or company purchasing the Products.

GOVERNING TERMS. Any purchase and shipment of Products shall be deemed to be on the terms and conditions stated herein. Any and all terms and conditions submitted by Purchaser are hereby rejected.

TAXES. The prices set forth herein do not include any sales, use, excise, ad valorem, property or other taxes applicable to the sale, use or delivery of the Products, all of which shall be paid by Purchaser separately or added to the contract price and paid by Purchaser to Keystone.

PRICES. The prices set forth for the Products are in United States Dollars. All prices shall be as of the invoice date, regardless of the date the Product is originally ordered from Keystone.

DELIVERY. Products sold hereunder shall be delivered F.O.B. shipping point. Purchaser shall assume responsibility for all subsequent delivery or shipping charges.

ACCEPTANCE. Purchaser must inspect all Products for defects and/or non-conformity within seven (7) calendar days of delivery. After this period, Purchaser will have accepted such Products.

PAYMENT. Payment terms shall be as set forth on the front of the invoice. If payment is late, the maximum statutory interest rate will apply on such late amount. Keystone may, in its sole discretion, suspend shipments of Products until full payment by Purchaser of all past due amounts to Keystone is made.

RETURN; EXCHANGE. Purchaser may return any Product within thirty (30) calendar days of the invoice date so long as there is an offsetting (exchange) order of equal or greater value. Purchaser may return any Product after thirty (30) calendar days from the invoice date, but no later than ninety (90) calendar days after the invoice date, so long as there is an offsetting (exchange) order of equal or greater value, provided, however, such returns shall be subject to a 15% re-stocking fee. Keystone will not accept any returns or exchanges after ninety (90) calendar days following the invoice date.

No biomaterial products (bone grafting and membrane) will be accepted for return or exchange. No returns or exchanges will be accepted for Products that (i) were originally sold as part of a kit, pack or set, (ii) were included on the original invoice at no charge, (iii) have less than six (6) months remaining before expiration or (iv) were custom ordered. The Purchaser shall pay all shipping costs associated with the return or exchange of a Product. No Products will be accepted by Keystone for return or exchange without an RMA number given to Purchaser by Keystone Customer Service. Returned Products must be unexpired, unopened and undamaged, returned freight prepaid.

WARRANTY. EXCEPT AS EXPRESSLY SET FORTH IN KEYSTONE'S LIFETIME IMPLANT WARRANTY PROGRAM, KEYSTONE DOES NOT MAKE ANY WARRANTY, EXPRESS OR IMPLIED, EXCEPT THAT ITS PRODUCTS SHALL BE FREE FROM DEFECTS IN MATERIAL AND/OR WORKMANSHIP. THIS WARRANTY APPLIES ONLY TO THE ORIGINAL PURCHASER. PURCHASER ASSUMES ALL RISKS AND LIABILITY RESULTING FROM THE USE OF THESE PRODUCTS, WHETHER USED SEPARATELY OR IN COMBINATION WITH OTHER PRODUCTS. THE WARRANTY DOES NOT APPLY IF THE FAILURE OF THE PRODUCT IS CAUSED BY TRAUMA OR ACCIDENT OR BY THE PATIENT IN WHOM THE IMPLANT WAS INSTALLED OR IF THE PURCHASER IS NOT CURRENT ON ALL OUTSTANDING PAYMENTS TO KEYSTONE AT THE TIME OF SUCH CLAIM. KEYSTONE MAY MODIFY OR TERMINATE ITS WARRANTY PROGRAM AT ANY TIME, HOWEVER ANY CHANGES TO THE WARRANTY PROGRAM WILL NOT AFFECT PRODUCTS ALREADY PURCHASED AT THE TIME OF SUCH CHANGE. PURCHASER'S SOLE REMEDY WITH RESPECT TO CLAIMS RELATED TO KEYSTONE PRODUCTS IS RESTRICTED TO THE WARRANTY POLICY.

LIMITATION OF LIABILITY. UNDER NO CIRCUMSTANCES SHALL KEYSTONE OR ANY OF ITS AFFILIATES BE LIABLE FOR BUSINESS INTERRUPTION, LOSS OF PROFITS, SPECIAL, INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE AND FROM ANY CAUSE WHETHER BASED IN CONTRACT, WARRANTY OR TORT (INCLUDING NEGLIGENCE) OR OTHER LEGAL THEORY, EVEN IF KEYSTONE HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

INTELLECTUAL PROPERTY. All intellectual property rights in Products sold to Purchaser shall remain the exclusive property of Keystone and its licensors.

PERSONAL DATA. It is Purchaser's responsibility to obtain patient's consent to transfer any personal data to Keystone.

TERMINATION. Keystone may terminate this or any agreement with Purchaser to purchase its products upon immediate written notice if Purchaser fails to pay any of its debt obligations to Keystone on time.

FORCE MAJEURE. Keystone shall not be liable for any loss or damage due to failure or delay arising out of any cause beyond the reasonable control of Keystone. In the event of any failure or delay resulting from such causes, an equitable adjustment of delivery and any other appropriate terms and conditions shall be made. No such failure or delay shall be the basis for an increase in Keystone's obligations, nor any termination by Purchaser.

VALIDITY. If any provision of these Terms and Conditions is found to be illegal or unenforceable in any respect, such illegality or unenforceability shall not affect any other provision of these Terms and Conditions, all of which shall remain enforceable in accordance with their terms.

GENERAL PROVISIONS. The purchase of the Products and these Terms and Conditions are governed by the laws of the Commonwealth of Massachusetts, without regard to conflicts of law principles. Any dispute, claim or controversy arising out of or in connection with the purchase of the Products shall be resolved by arbitration in Boston, Massachusetts. This document incorporates all oral and written representations between the parties and constitutes the entire agreement and understanding of the parties with respect to the subject matter hereof and supersedes any and all other agreements either oral or written between the parties with respect to such subject matter. No amendment or modification may be made to this document unless in writing and duly executed by an authorized representative of Keystone. Keystone may subcontract, assign or transfer its obligations or rights to a competent third party. Purchaser may not assign or transfer any of its obligations without the consent of Keystone.

Keystone Dental, Inc.

154 Middlesex Turnpike
Burlington, MA 01803 USA
General Inquiries: 781-328-3600
Customer Relations: 866-902-9272 or 781-328-3490
Fax: 866-903-9272 or 781-328-3400
Email: info@keystonedental.com



**Keystone Dental SpA
European Headquarters**

via A. Fleming, 19
37135 Verona, Italy
Phone: +39 045 8230294
Fax: +39 045 8250296
Email: info.it@keystonedental.com



Caution, consult accompanying documents

Trademark Acknowledgements

TiLobe, Quick-Cap and "Smarter Thinking. Simpler Design."
are trademarks and registered trademarks of Keystone Dental, Inc.